



## Universal Fit Adjustable Height Rear Anti-Tipper

This innovative design eliminates piles of mismatched anti-tippers because these fit 3/4", 7/8" or 1" tubing. They feature infinite height adjustment and 2" easy gliding wheels. No more too high, too low or wrong mounting style!

Universal Anti-Tipper    193-151    \$42.86 a pair

## Universal Fit Adjustable Height Rear Anti-Tipper

Made of strong bar stock, this anti-tipper will help you streamline your inventory. It's the only anti-tipper that you need to stock. New quick adjust feature allows the anti-tipper to be positioned properly without the use of tools. Flat black powder coat finish or Chrome finish available. 2 Inch Durable Plastic Wheel.

Black Anti-Tipper    193-171    \$34.82 a pair  
 Chrome Anti-Tipper    193-172    \$34.82 a pair



## Adjustable Anti-Theft Bar

This device prevents the wheelchair from folding up. It works on wheelchairs from 16" to 24" wide.



### Anti-Theft Bar

16" to 20" Wide    199-109    \$15.08  
 22" to 24" Wide    199-110    \$16.47

## Over Head Anti-Theft Bar

This over head device prevents the wheelchair from folding up. Fits wheelchairs 18" to 26" wide. Other sizes available upon request. No available in non-magnetic for MRI wheelchairs.



### STD Wheelchairs

18" Wide    199-111    \$50.00 ea.  
 20" Wide    199-115    \$83.40 ea.  
 22" Wide    199-118    \$83.36 ea.  
 24" Wide    199-116    \$83.36 ea.  
 26" Wide    199-117    \$83.36 ea.

### MRI Wheelchair (non-magnetic)

18" Wide    199-130    \$135.27 ea.  
 20" Wide    199-132    \$171.68 ea.  
 22" Wide    199-131    \$171.68 ea.  
 24" Wide    199-133    \$171.68 ea.

## Universal Wheelchair Brake Extensions

### Features:

- > Available in three lengths
- > Available in two slot sizes
- > Available with or without elastic cord

### Elastic Cord

The elastic cord connects the extension to the brake. This helps to hold the extension on the brake.

- A Without Elastic Cord**    \$4.95 ea.
- 6" Extension with 1/8" Slot Size    154-102
  - 6" Extension with 3/16" Slot Size    154-105
  - 8" Extension with 1/8" Slot Size    154-101
  - 8" Extension with 3/16" Slot Size    154-104
  - 10" Extension with 1/8" Slot Size    154-100
  - 10" Extension with 3/16" Slot Size    154-103
- B With Elastic Cord**    \$7.95 ea.
- 6" Extension with 1/8" Slot Size    154-108
  - 6" Extension with 3/16" Slot Size    154-111
  - 8" Extension with 1/8" Slot Size    154-107
  - 8" Extension with 3/16" Slot Size    154-110
  - 10" Extension with 1/8" Slot Size    154-106
  - 10" Extension with 3/16" Slot Size    154-109

