

Grey Mag Heavy Duty



24" x 1" 8 Spoke Grey Mag, Solid Grey Urethane Tire with $\frac{3}{4}$ " Chrome Handrim. Hub width measured Outside of Bearing to Outside of Bearing is approximately $2\frac{5}{8}$ ". Has flush outer bearing.

Bearing I.D. $\frac{5}{8}$ " 160-415 \$75.50 ea. Everest & Jennings

Pneumatic Tire



24" x 1 $\frac{1}{4}$ " 8 Spoke Black or Grey Mag, Grey Urethane Full Profile Treaded Tire with $\frac{3}{4}$ " Chrome Handrim. Hub width measured Outside of Bearing to Outside of Bearing is approximately $2\frac{5}{8}$ ". Has flush outer bearing.

Bearing I.D. $\frac{7}{16}$ " 160-607 \$80.00 ea. Black Mag
 $\frac{5}{8}$ " 160-608 \$80.00 ea. Black Mag
 $\frac{1}{2}$ " 160-609 \$80.00 ea. Black Mag

Bearing I.D. $\frac{7}{16}$ " 160-407 Disc. Grey Mag
 $\frac{5}{8}$ " 160-408 Disc. Grey Mag
 $\frac{1}{2}$ " 160-409 Disc. Grey Mag

Everest & Jennings
Theradyne, Quickie,
Invacare & Others



24" x 1 $\frac{3}{8}$ " 8 Spoke Grey Mag, Grey Pneumatic Tire with $\frac{3}{4}$ " Chrome Handrim. Hub width measured Outside of Bearing to Outside of Bearing is approximately $2\frac{5}{8}$ ". Has flush outer bearing.

Bearing I.D. $\frac{7}{16}$ " 160-911 \$54.95 ea.
 $\frac{5}{8}$ " 160-912 \$54.95 ea.
 $\frac{1}{2}$ " 160-913 \$54.95 ea.

Everest & Jennings
Theradyne, Quickie,
Invacare & Others



24" x 1 $\frac{3}{8}$ " 8 Spoke Black Mag, Grey Pneumatic Tire with $\frac{3}{4}$ " Chrome Handrim. Hub width measured Outside of Bearing to Outside of Bearing is approximately $2\frac{5}{8}$ ". Has flush outer bearing.

Bearing I.D. $\frac{7}{16}$ " 160-611 \$60.85 ea.
 $\frac{5}{8}$ " 160-612 \$60.85 ea.
 $\frac{1}{2}$ " 160-613 \$60.85 ea.

Everest & Jennings
Theradyne, Quickie,
Invacare & Others



24" x 1 $\frac{3}{4}$ " - 2 $\frac{1}{8}$ " 7 Spoke Black Mag, Grey Pneumatic Tire with Stainless Steel Handrim. Hub width measured Outside of Bearing to Outside of Bearing is approximately $1\frac{7}{16}$ ". Has recessed outer bearing.

Bearing I.D. $\frac{1}{2}$ " with $\frac{7}{16}$ " Reducer 168-711 \$98.35 ea.

Theradyne, Quickie, Invacare,
Quadra & Others

MRI Non-Magnetic Rear Wheel, Complete Assembly



Rear Wheel for 18" to 20" Wide Standard MRI Wheelchair with $\frac{7}{16}$ " Bearings
Includes Wheel, Bearings, Axle and Nut

160-950 \$216.16 ea.



Rear Wheel for 22" to 24" Wide Standard MRI Wheelchair with $\frac{5}{8}$ " Bearings
Includes Wheel, Bearings, Axle and Nut

160-951 \$216.16 ea.



Rear Wheel for 24" to 26" Wide Heavy Duty MRI Wheelchair with $\frac{5}{8}$ " Bearings
Includes Wheel, Bearings, Axle and Nut

160-952 \$232.42 ea.

One Arm Drives

For ordering one arm drives call toll free **800-328-5343** with the following information:

Make of Chair
Model of Chair
Serial Number

