

## Wheelchair I.V. Poles and Oxygen Holders

**Now Available for MRI Wheelchairs  
Non-Magnetic for Use In MRI Room**

### A I.V. Pole, 2 Hook

*Fits most standard chairs, I.V. pole is adjustable*

STD Chair 198-100 \$31.00

MRI Chair (non-magnetic) 198-110 \$123.67

### B I.V. Pole, 1 Hook

*Fits most standard chairs, I.V. pole is adjustable*

STD Chair 198-101 \$67.50

### C Oxygen Tank Holder

*Fits most standard chairs, oxygen holder is for D-E size cylinders*

STD Chair Right Side Mount 198-200 \$48.39

STD Chair Left Side Mount 198-201 \$48.39

MRI Chair (non-magnetic) 198-210 \$109.83

Right Side Mount

### D I.V. Pole and Oxygen Tank Holder

*Fits most standard chairs, adjustable in height, oxygen holder is for D-E size cylinders*


STD Chair 198-300 \$59.95

MRI Chair (non-magnetic) 198-310 \$182.30

### E O2 Tank Carrier for Rollators

*Attaches conveniently to the frame of a 3- or 4-wheel rollator. Provides secure handling of a single D or E cylinder where a mobile oxygen unit is required. Fits most rollators*

STD Rollators 250-193 \$26.67



## Front Anti-Tippers

This product prevents patients from tipping forward on the wheelchair. Fits on 7/8" and 1" O.D. tubing.

### Front Anti-Tippers

7/8" O.D. Tubing 193-102 \$44.89

1" O.D. Tubing 193-119



## Adjustable Front Anti-tippers

### Features:

- > Fits both 7/8" and 1" frame diameter
- > Made from durable powder coated steel tubing and heavy gauge chrome plated steel
- > 2" easy gliding wheels
- > Mounting hardware and easy to follow installation instructions included
- > Five height adjustments; fits full height to low hemi chairs

