

Wheels

Toll Free 800-328-5343

Fax Free 866-829-9744

www.ocelco.com

Our wheels are available in many configurations including Grey and black mags, pneumatic and solid black or Grey tires, varying hub widths (be sure to consider your axle length when determining the hub size you need), multiple numbers of spokes, many different handrim styles and different bearing inside diameters.



20" x 1³/₈" Black Mag, Pneumatic Tire with ³/₄" Powder Coated Handrim. Hub width measured Outside of Bearing to Outside of Bearing is approximately 2³/₈". Has flush outer bearing.

Bearing I.D. 7/16" 168-211 \$81.62 ea.

**Everest & Jennings,
Theradyne, Quickie,
Invacare & Others**



20" x 1³/₈" Black Mag, Pneumatic Tire with ³/₄" Stainless Steel Handrim. Hub width measured Outside of Bearing to Outside of Bearing is approximately 2³/₈". Has recessed outer bearing.

Bearing I.D. 1/2" with 7/16" Reducer 168-311 \$85.03 ea.

Invacare & Quickie



20" x 1³/₄" 5 Spoke Black Mag, Full Profile Urethane Tire with ³/₄" Aluminum Powder Coated Handrim. Hub width measured Outside of Bearing to Outside of Bearing is approximately 2⁵/₈". Has flush outer bearing.

Bearing I.D. 7/16" 168-411 \$86.34 ea.
 5/8" 168-412 \$86.34 ea.

**Everest & Jennings,
Theradyne, Quickie,
Invacare & Others**



22" x 1" 8 Spoke Black Mag, Low Profile Black Urethane Tire with ³/₄" Aluminum Powder Coated Handrim. Hub width measured Outside of Bearing to Outside of Bearing is approximately 2⁵/₈". Has flush outer bearing.

Bearing I.D. 7/16" 168-011 \$84.17 ea.
 5/8" 168-012 \$82.03 ea.
 1/2" 168-013 \$90.13 ea.

**Everest & Jennings
Theradyne, Quickie,
Invacare & Others**



22" x 1³/₈" 7 Spoke Black Mag, Pneumatic Tire with ³/₄" Handrim. Hub width measured Outside of Bearing to Outside of Bearing is approximately 2⁵/₈". Has flush outer bearing.

Bearing I.D. 7/16" 168-111 \$82.97 ea.
 5/8" 168-112 \$82.31 ea.
 1/2" 168-113 \$90.40 ea.

**Everest & Jennings
Theradyne, Quickie,
Invacare & Others**



22" x 1³/₈" 8 Spoke Black Mag, Full Profile Poly Tire with ³/₄" Aluminum Powder Coated Handrim. Has flush outer bearing and fits most **Invacare, E&J, Theradyne and Quickie** chairs.

Bearing I.D. 7/16" 168-511 \$85.61 ea. Hub width measured Outside of Bearing to Outside of Bearing is approximately 2⁵/₈"
 1/2" 168-512 \$89.58 ea. Hub width measured Outside of Bearing to Outside of Bearing is approximately 2¹/₁₆"

**Everest & Jennings
Theradyne, Quickie,
Invacare & Others**



24" x 1" Alloy Spoke Rim, Solid Grey Tire with ⁵/₈" Chrome Handrim. Has flush outer bearing. 2¹/₂" Hub width.

Bearing I.D. 7/16" 160-111 \$40.83 ea.
 5/8" 160-112 \$44.61 ea.
 1/2" 160-113 \$49.95 ea.

**Everest & Jennings
Theradyne, Quickie,
Invacare & Others**

Black Mag



24" x 1" 8 Spoke Black Mag, Solid Grey Urethane Tire with ³/₄" Handrim. Hub width measured Outside of Bearing to Outside of Bearing is approximately 2⁵/₈". Has flush outer bearing.

Bearing I.D.	7/16"	160-411	\$65.74 ea.	Everest & Jennings Theradyne, Quickie & Others
	5/8"	160-412	\$65.74 ea.	
	1/2"	160-413	\$65.74 ea.	
	12mm x 32mm	160-414	\$65.74 ea.	



24" x 1" 8 Spoke Black Mag, Black Urethane Low Profile Tire with ³/₄" Chrome Handrim. Hub width measured Outside of Bearing to Outside of Bearing is approximately 2¹/₄". Has recessed outer bearing.

Bearing I.D.	7/16"	160-507	\$81.34 ea.	Invacare
	5/8"	160-508	\$81.34 ea.	
	1/2"	160-509	\$81.34 ea.	



24" x 1" 8 Spoke Black Mag, Solid Black Urethane Tire with ³/₄" Chrome Handrim. Hub width measured Outside of Bearing to Outside of Bearing is approximately 2⁵/₈". Has flush outer bearing.

Bearing I.D.	7/16"	160-811	\$57.95 ea.	Everest & Jennings Theradyne, Quickie, Invacare & Others
	5/8"	160-812	\$57.95 ea.	
	1/2"	160-813	\$57.95 ea.	



24" x 1" 9 Spoke Black Mag, Solid Urethane Tire with ³/₄" Chrome Handrim. Hub width measured Outside of Bearing to Outside of Bearing is approximately 2¹/₄". Has recessed outer bearing.

Bearing I.D.	7/16"	160-711	\$49.95 ea.	Grey DTM Tire	Older Invacare Chairs
	5/8"	160-712	\$49.95 ea.	Grey DTM Tire	
	1/2"	160-713	\$49.95 ea.	Grey DTM Tire	



Bearing I.D.	7/16"	168-611	\$70.57 ea.	Black Tire	New Invacare Chairs
	5/8"	168-612	\$70.72 ea.	Black Tire	
	1/2"	168-613	\$70.72 ea.	Black Tire	
Drive Rear Wheel Assembly, Black Handrim, Hub Width 2¹/₄".					Drive Chairs, Cruiser and Silver Sport
35 mm	160-114	\$33.76	Black Tire		

Grey Mag



24" x 1" 7 Spoke Grey Mag, Solid Grey Urethane Tire with ³/₄" Grey Plastic Handrim. Hub width measured Outside of Bearing to Outside of Bearing is approximately 2⁵/₈". Has flush outer bearing.

Bearing I.D.	7/16"	160-511	\$57.95 ea.	Everest & Jennings Theradyne, Quickie & Others
	5/8"	160-512	\$57.95 ea.	
	1/2"	160-513	\$57.95 ea.	



24" x 1" 9 Spoke Grey Mag, Solid Urethane Tire with ³/₄" Chrome Handrim. Hub width measured Outside of Bearing to Outside of Bearing is approximately 2¹/₄". Has recessed outer bearing.

Bearing I.D.	7/16"	160-011	\$49.95 ea.	Grey Tire	Invacare
	5/8"	160-012	\$49.95 ea.	Grey Tire	
	1/2"	160-013	\$49.95 ea.	Grey Tire	



24" x 1" 8 Spoke Grey Mag, Solid Dark Grey Urethane Tire with ³/₄" Chrome Handrim. Hub width measured Outside of Bearing to Outside of Bearing is approximately 2⁵/₈". Has flush outer bearing.

Bearing I.D.	7/16"	160-211	\$57.59 ea.	Everest & Jennings Theradyne, Quickie, Invacare & Others
	5/8"	160-212	\$57.59 ea.	
	1/2"	160-213	\$57.59 ea.	



Grey Mag Heavy Duty

24" x 1" 8 Spoke Grey Mag, Solid Grey Urethane Tire with $\frac{3}{4}$ " Chrome Handrim. Hub width measured Outside of Bearing to Outside of Bearing is approximately $2\frac{5}{8}$ ". Has flush outer bearing.

Bearing I.D. $\frac{5}{8}$ " 160-415 \$75.50 ea. Everest & Jennings



Pneumatic Tire

24" x 1 $\frac{1}{4}$ " 8 Spoke Black or Grey Mag, Grey Urethane Full Profile Treaded Tire with $\frac{3}{4}$ " Chrome Handrim. Hub width measured Outside of Bearing to Outside of Bearing is approximately $2\frac{5}{8}$ ". Has flush outer bearing.

Bearing I.D. $\frac{7}{16}$ " 160-607 \$80.00 ea. Black Mag
 $\frac{5}{8}$ " 160-608 \$80.00 ea. Black Mag
 $\frac{1}{2}$ " 160-609 \$80.00 ea. Black Mag

Bearing I.D. $\frac{7}{16}$ " 160-407 Disc. Grey Mag
 $\frac{5}{8}$ " 160-408 Disc. Grey Mag
 $\frac{1}{2}$ " 160-409 Disc. Grey Mag

Everest & Jennings
Theradyne, Quickie,
Invacare & Others



24" x 1 $\frac{3}{8}$ " 8 Spoke Grey Mag, Grey Pneumatic Tire with $\frac{3}{4}$ " Chrome Handrim. Hub width measured Outside of Bearing to Outside of Bearing is approximately $2\frac{5}{8}$ ". Has flush outer bearing.

Bearing I.D. $\frac{7}{16}$ " 160-911 \$54.95 ea.
 $\frac{5}{8}$ " 160-912 \$54.95 ea.
 $\frac{1}{2}$ " 160-913 \$54.95 ea.

Everest & Jennings
Theradyne, Quickie,
Invacare & Others



24" x 1 $\frac{3}{8}$ " 8 Spoke Black Mag, Grey Pneumatic Tire with $\frac{3}{4}$ " Chrome Handrim. Hub width measured Outside of Bearing to Outside of Bearing is approximately $2\frac{5}{8}$ ". Has flush outer bearing.

Bearing I.D. $\frac{7}{16}$ " 160-611 \$60.85 ea.
 $\frac{5}{8}$ " 160-612 \$60.85 ea.
 $\frac{1}{2}$ " 160-613 \$60.85 ea.

Everest & Jennings
Theradyne, Quickie,
Invacare & Others



24" x 1 $\frac{3}{4}$ " - 2 $\frac{1}{8}$ " 7 Spoke Black Mag, Grey Pneumatic Tire with Stainless Steel Handrim. Hub width measured Outside of Bearing to Outside of Bearing is approximately $1\frac{7}{16}$ ". Has recessed outer bearing.

Bearing I.D. $\frac{1}{2}$ " with $\frac{7}{16}$ " Reducer 168-711 \$98.35 ea.

Theradyne, Quickie, Invacare,
Quadna & Others



MRI Non-Magnetic Rear Wheel, Complete Assembly

Rear Wheel for 18" to 20" Wide Standard MRI Wheelchair with $\frac{7}{16}$ " Bearings
Includes Wheel, Bearings, Axle and Nut
160-950 \$216.16 ea.

Rear Wheel for 22" to 24" Wide Standard MRI Wheelchair with $\frac{5}{8}$ " Bearings
Includes Wheel, Bearings, Axle and Nut
160-951 \$216.16 ea.


Rear Wheel for 24" to 26" Wide Heavy Duty MRI Wheelchair with $\frac{5}{8}$ " Bearings
Includes Wheel, Bearings, Axle and Nut
160-952 \$232.42 ea.



One Arm Drives

For ordering one arm drives call toll free **800-328-5343** with the following information:

Make of Chair
Model of Chair
Serial Number



All Hub Widths are approximate and measure outside of bearing to outside of bearing.



5" x 1"

3 Spoke Black Mag, Black Urethane Tire.

Bearing I.D.	Hub Width	Part No.	Price
5/16" Precision Bearing	1"	161-600	\$25.06 ea.

Solid Tire



8 Spoke Black Mag, Black Urethane Tire.

Bearing I.D.	Hub Width	Part No.	Price
5/16" Precision Bearing	1 1/2"	161-615	\$23.95 ea.
7/16" Flange Bearing	1 3/4"	161-608	\$23.95 ea.
5/16" Flange Bearing	1 3/4"	161-601	\$23.95 ea.

Solid Tire



6" x 1"

8 Spoke Black Mag, Black Rubber Tire.

Bearing I.D.	Hub Width	Part No.	Price
5/16" Flange Bearing	1 1/2"	161-452	\$30.34 ea.

Solid Tire



6 Spoke Black Mag, Grey Urethane Tire.

Bearing I.D.	Hub Width	Part No.	Price
5/16" Precision Bearing	1"	161-616	\$26.97 ea.

Solid Tire



8 Spoke Black Mag, Black Urethane Tire.

Bearing I.D.	Hub Width	Part No.	Price
5/16" Precision Bearing	1 1/2"	161-617	\$25.40 ea.
7/16" Flange Bearing	1 3/4"	161-609	\$25.40 ea.
5/16" Flange Bearing	1 3/4"	161-602	\$25.40 ea.

Solid Tire



6" x 1 1/4"

3 Spoke Black Mag, Grey Pneumatic Tire.

Bearing I.D.	Hub Width	Part No.	Price
5/16" Precision Bearing	1"	161-462	\$40.80 ea.
7/16" Flange Bearing	1 1/4"	161-483	\$33.93 ea.

Pneumatic Tire



8 Spoke Black Mag, Black Urethane Tire.

Bearing I.D.	Hub Width	Part No.	Price
5/16" Precision Bearing	1 1/2"	161-618	\$23.95 ea.
7/16" Flange Bearing	1 1/2"	161-610	\$23.95 ea.
5/16" Flange Bearing	1 1/2"	161-603	\$23.95 ea.

Solid Tire



7" x 1"

8 Spoke Black Mag, Black Urethane Tire.

Bearing I.D.	Hub Width	Part No.	Price
5/16" Precision Bearing	1 1/2"	161-619	\$32.80 ea.
7/16" Flange Bearing	1 1/2"	161-611	\$36.35 ea.
5/16" Flange Bearing	1 1/2"	161-604	\$35.61 ea.

Solid Tire



7 1/2" x 1"

8 Spoke Black Mag, Black Urethane Tire.

Bearing I.D.	Hub Width	Part No.	Price
5/16" Precision Bearing	1 1/2"	161-620	\$34.44 ea.
7/16" Flange Bearing	2 3/8"	161-612	\$25.13 ea.
5/16" Flange Bearing	1 1/2"	161-605	Disc.

Solid Tire
Black Mag Shown



8" x 1"

7 Spoke, Grey Rubber Tire.

Bearing I.D.	Hub Width	Part No.	Price
*5/16" Bearing Black Mag	2 5/8"	161-412	\$15.95 ea.
*5/16" Bearing Grey Mag	2 5/8"	161-212	\$15.95 ea.

*Standard Everest & Jennings

Solid Tire
Black Mag Shown



8 Spoke, Grey Rubber Tire.

Bearing I.D.	Hub Width	Part No.	Price
7/16" Bearing Grey Mag	2 3/8"	161-211	\$15.95 ea.
*7/16" Bearing Black Mag	2 3/8"	161-411	\$15.95 ea.
1/4" Bearing Black Mag	2 3/8"	161-413	\$19.95 ea.
1/4" Bearing Grey Mag	2 3/8"	161-213	\$19.95 ea.

Solid Tire
Black Mag Shown

*Old Standard Invacare



5 Spoke Heavy Duty Mag, Grey Rubber Tire.

Bearing I.D.	Hub Width	Part No.	Price
5/16" Bearing Grey Mag	2 5/16"	161-312	\$23.55 ea.
5/16" Bearing Black Mag	2 1/2"	161-512	\$27.95 ea.
7/16" Bearing Grey Mag	2 3/8"	161-311	\$23.55 ea.
7/16" Bearing Black Mag	2 3/8"	161-511	\$27.95 ea.

Solid Tire



8 Spoke Black Mag, Black Urethane Tire.

Bearing I.D.	Hub Width	Part No.	Price
5/16" Precision Bearing	1 1/2"	161-621	\$23.95 ea.
*7/16" Flange Bearing	1 3/4"	161-613	\$23.95 ea.
5/16" Flange Bearing	1 3/4"	161-606	\$23.95 ea.
8 mm Flange Bearing	56 mm	161-414	\$13.55 ea.

Solid Tire



Wire Spoke Wheel, Grey Rubber Tire.

Bearing I.D.	Hub Width	Part No.	Price
1/2" Flange Bearing	2 3/8"	161-114	\$34.95 ea.

Solid Tire



8" x 1"

MRI Non-Magnetic

Front Wheel for 18" to 24" Wide Standard Wheelchair

161-640 \$30.17 ea.

8" x 1" Wheel Assembly for 18" to 24" MRI Chairs

Includes fork, wheel, bearings, axle, and nut. 161-641 \$94.21 ea.

All Hub Widths are approximate and measure outside of bearing to outside of bearing.

8" x 1 1/4"



8 Spoke Black Mag, Grey Pneumatic Tire.

Bearing I.D.	Hub Width		
5/16" Bearing	1 3/4"	161-472	Disc.
7/16" Bearing	2 3/8"	161-471	\$26.78 ea.

Pneumatic Tire



5 Spoke Grey Mag, Grey Pneumatic Tire.

Bearing I.D.	Hub Width		
5/16" Bearing	2 3/8"	161-262	\$39.75 ea.
7/16" Bearing	2 1/8"	161-261	Disc.

Pneumatic Tire



7 Spoke Black Mag, Black Urethane Tire.

Bearing I.D.	Hub Width		
5/16" Bearing	1 3/4"	161-491	\$45.00 ea.

Solid Tire



8 Spoke Black Mag, Grey Pneumatic Tire with Flat Free Poly-Foam Insert.

Bearing I.D.	Hub Width		
5/16" Bearing	1 3/4"	161-493	\$50.53 ea.
7/16" Bearing	1 3/4"	161-494	\$50.53 ea.

Pneumatic Flat Free

8" x 1 3/4"



2 Piece Black Hub, Grey Pneumatic Tire.

Bearing I.D.	Hub Width		
5/16" Bearing	2 1/2"	161-422	Disc.
7/16" Bearing	2 1/2"	161-421	Disc.

Pneumatic Tire



2 Piece Black Hub, Black Urethane Tire.

Bearing I.D.	Hub Width		
5/16" Bearing	2 1/2"	161-432	Disc.
7/16" Bearing	2 1/2"	161-431	Disc.

Solid Tire



2 Piece Grey Hub, Grey Pneumatic Tire.

Bearing I.D.	Hub Width		
5/16" Bearing	2 1/2"	161-222	\$62.90 ea.
7/16" Bearing	2 1/2"	161-221	\$57.60 ea.

Pneumatic Tire



8" x 1 3/4"
2 Piece Grey Hub, Dark Grey Urethane Tire.

Bearing I.D.	Hub Width		
5/16" Bearing	2 1/2"	161-232	Disc.
7/16" Bearing	2 1/2"	161-231	Disc.

Solid Tire



2 Piece Aluminum Hub, Grey Pneumatic Tire with Flat Free Poly-Foam Insert.

Bearing I.D.	Hub Width		
5/16" Bearing	2 1/2"	161-495	Disc.
7/16" Bearing	2 1/2"	161-496	\$51.16 ea.

Pneumatic Flat Free

8" x 2"



3 Spoke Grey Mag, Grey Pneumatic Tire.

Bearing I.D.	Hub Width		
5/16" Bearing	2 1/2"	161-223	\$27.10 ea.
7/16" Bearing	2 1/2"	161-224	\$27.10 ea.

Pneumatic Tire



2 Piece Black Hub, Grey Pneumatic Tire.

Bearing I.D.	Hub Width		
5/16" Bearing	2 1/2"	161-482	\$27.93 ea.
7/16" Bearing	2 1/2"	161-481	\$28.67 ea.

Pneumatic Tire



2 Piece Black Hub, Black Urethane Tire.

Bearing I.D.	Hub Width		
5/16" Bearing	2 1/2"	161-607	\$28.10 ea.
7/16" Bearing	2 1/2"	161-614	\$28.10 ea.

Solid Tire



2 Piece Black Hub, Grey Flat Free Pneumatic Poly Foam Filled Tire.

Bearing I.D.	Hub Width		
5/16" Bearing	2 1/2"	161-497	\$44.17 ea.
7/16" Bearing	2 1/2"	161-498	\$43.19 ea.

Pneumatic Flat Free

8" x 2"

MRI Non-Magnetic



Front Wheel for 24" to 26" Wide Heavy Duty MRI Wheelchair
 161-642 \$47.78 ea.



8" x 2" Wheel Assembly for 24" to 26" MRI Heavy Duty Chairs
 Includes fork, wheel, bearings, axle, and nut. 161-643 \$110.43 ea.