



New style safety buckle



Safety Belt

Made with long wearing web materials. Choose between a simple two piece belt with velcro buckle, plastic buckle or auto buckle. Equipped with your choice of ends to attach to wheelchair, axle mount or beneath upholstery screws. The belts come in three standard sizes, special sizes are available upon request. Specify color of webbing when ordering your belts.

Note: Length of belt includes 18" for back of wheelchair. Example: Actual length of 60" belt is only 42".

Velcro Buckle

A. Upholstery Screw Mount	\$12.05	B. Axle Mount	\$12.85
60" Belt	190-111-CL	60" Belt	190-129-CL
72" Belt	190-103-CL	72" Belt	190-130-CL
90" Belt	(\$12.75) 190-113-CL	90" Belt	190-131-CL

Plastic Buckle

C. Upholstery Screw Mount	\$12.55	D. Axle Mount	\$13.35
60" Belt	190-112-CL	60" Belt	190-126-CL
72" Belt	190-124-CL	72" Belt	190-127-CL
90" Belt	190-125-CL	90" Belt	190-128-CL

Auto Buckle

E. Upholstery Screw Mount	\$19.10	F. Axle Mount	\$19.90
60" Belt	190-104-CL	60" Belt	190-116-CL
72" Belt	190-122-CL	72" Belt	190-120-CL
90" Belt	190-123-CL	90" Belt	190-121-CL

See Page 159 for Web Colors



Safety Belt & Cuff

The patented cuff provides added safety. Eliminates easy access to the velcro closure ensuring fewer accidents. Mounts beneath upholstery screws on back of chair. Special sizes upon request. See Page 159 for Web Colors

Safety Belts and Cuff	\$12.10
Wheelchair Belt 60"	190-100-CL
Wheelchair Belt 72"	190-101-CL
Geri-Chair 72"	190-102-CL



Butterfly Safety Belt

Made with 1½" web, 1½" plastic side release buckles and Naugahyde. Attaches to wheelchair beneath upholstery screws or screws to solid backs. Ideal for patients needing extra support while sitting in chair. Specify color of Naugahyde and webbing when ordering. Special sizes available. See Page 159 for Naugahyde Colors

Butterfly Safety Belt	\$31.65
Small	190-117-CL
Medium	190-118-CL
Large	190-119-CL

Alarmed Velcro Seat Belt

This product is ideal for the patient who needs the reminder of a velcro lap belt but who, upon request, is able to release the belt. At the instant the belt is separated, an alarm sounds reminding the patient to remain seated. The sound of the alarm also alerts the staff of a potential fall. The seat belt features a red pull strap to facilitate locating and releasing the belt. The rugged alarm unit features 95db. horn and stainless steel mounting clip. A 9 volt battery is included.



Seat Belt Medium 36" to 46" 193-134 \$85.71

Seat Belt Large 46" to 58" 193-134-LG \$85.71

Seat Belt X-Large 58" to 70" 193-134-XL \$85.71

Fall Prevention Monitor

The Safe-t Mate fall prevention monitor has been designed to overcome the pitfall of others on the market. Many monitors are expensive yet unreliable. One fall and they break. Some are less costly but are not designed or intended to be used as patient fall monitors. This monitor incorporates durability & versatility into an attractive, remarkably economical package. Features include: 1) Custom stainless clip that mounts to bed, wheelchair or walker without the need for additional mounting hardware, 2) Slide cord that adjusts for chair or bed use, 3) Pull away magnetic disc design allows for a variety of mounting positions. A 9 volt battery is included.



Fall Prevention Monitor 193-132 \$25.00

Spare Parts Kit 193-138 Disc.

UMP Chair Sentry Standard and Deluxe Monitors

Features:

- › Dual button reset; must be pressed simultaneously to silence alarm
- › Volume and delay switches inside battery compartment
- › No on/off switch
- › Self-timing pad use with monitor enables visual and audible indicators for pad expiration
- › Built-in Sentry cord stain relief
- › Low battery indicator, Uses single 9V battery
- › Selectable 4-tone alarm
- › Deluxe model has an integrated cushion protector

Sentry Standard Monitors

Chair Monitor 305-314 \$51.68

Cushion Protector 305-351 \$13.25

Sentry Deluxe Monitors

Chair Monitor 305-311 \$57.42



Standard Monitor



Deluxe Monitor



Cushion Protector

UMP Standard Chair Pads

Features:

- Durable vinyl cover cleans easily
- Compatible with all UMP monitors
- Available in 45 day or 1 year

**Other pads
available online at
www.ocelco.com or
www.WheelchairParts.com**

Ump Standard Chair Pads

45 Day Pad 305-321-45 \$19.97

45 Day Pad 6 Pack 305-329-45 \$102.69

1 Year Pad 305-321-01 \$43.27

1 Year Pad 6 Pack 305-329-01 \$204.95

